

Erasmus MC

University Medical Center Rotterdam



Reverse phenotyping: MCD Detection without clinical diagnosis

Diagnostic: NGS panels versus WES

Martina Wilke

Laboratory Specialist Clinical Genetics


Dept. of Clinical Genetics

Erasmus MC Rotterdam

The Netherlands


Milan, 19 June 2018

Declaration of Conflict of Interest




“I have no commercial disclosure”

Malformations of Cortical Development (MCD)




<u>Processes</u>	<u>Abnormalities</u>	<u>Time</u>
1. Dorsal Induction	Neural Tube Defects	3-7 w.
2. Ventral Induction	Holoprosencefalie	5-6 w.
3. Neuronal Glial Proliferation	Micro/Megalencefalie	8-16 w.
4. Migration	Lissencefalie en heterotopia	12-20 w.
5. Organization	Polymicrogyrie, dysplasie	>24 w.
6. Myelinatie	Hypo/dysmyelinatie	>24w.-2 yr.

The past: Sanger-sequencing



- **(1) Abnormal proliferation of the cortex (micro- en macrocephaly)**
 - ERCC2, ERCC5, ERCC6, CEP152, EIF2AK3, WDR62, TSC1/TSC2, AKT3
 - **(2) Abnormal neuronal migration (lissencephalie/heterotopie)**
PAFAH1B1/LIS1 (+ **MLPA**), TUBA1A, DCX, RELN, ARX, WDR62, FLNA, ARFGEF2
 - **(3) Abnormal organisation of the cortex and cerebrovascular abnormalities**
 - KIAA1279 / KBP, RTTN, GPR56, COL4A1, COL4A2
 - **Disadvantage:** (1) Time consuming; (2) Expensive
- *Still important: clear defined phenotype*




Single gene approach: 1227 patients

Gene	Nr of patients
DCX	14/ 84
ARX	3/ 177
ARFGEF2	5/
PAFAH1B1	24
TUBA1A	6
FLNA	33
COL4A1	27/ 422
COL4A2	1
ADGRG1	3
KIF1BP	3
ERCC2	2
ERCC5	2
EIF2AK3	1

Total of **124** mutations in 13 different genes

> **200** VUS

Diagnostic yield: Single gene approach




1227 individuals with cortical brain malformations

Direct diagnosis: **89 (7,2%)**

After evaluation: **35 (2,8%)**

Total diagnostic yield: **10 %**

Next Generation sequencing



From 22 → 103 genes

(1) Abnormal proliferation of the cortex (micro- en macrocephaly)

ASPM, MCPH1, STIL, CDK5RAP2, CENPJ, ATR, CEP152, PCNT, ERCC1, ERCC2, ERCC5, ERCC6, SLC25A19, CASK, TSEN54, EMG1, RNU4ATAC / U4actc snRNA, NBS1, MYCN, JAM3, PNKP, MED17, VPS13B/COH1, DKC1, EIF2AK3, IER3IP1, ASXL1, LYK5/STRADA, TSC1/ TSC2, PTEN, AKT1, AKT3, PIK3R2, PIK3CA, TBC1D7, TUBB5, KIF5C, STAMBP, TBC1D20

(2) Abnormal neuronal migration (lissencephalie/heterotopia)


PAFAH1B1/LIS1, TUBA1A, DCX, RELN, ARX, NDE1, WDR62, ACTB, ACTG1, LAMC3, FLNA, ARFGEF2, INTS8, YWHAE, ERMARD, DCHS1, FAT4, KIF2A, TUBG1

(3) Abnormal organisation of the cortex and cerebrovascular abnormalities

TUBB2B, TUBB3, TUBA8, EOMES/TBR2, TBC1D24, OCLN, RAB3GAP1, RAB3GAP2, RAB18, CNTNAP2, SNAP29, VLDLR, SRPX2, COL18A1, PAX6, KIAA1279/KBP, ATP6V0A2, RTTN, GPR56, FKRP, LARGE, POMGnT1, POMT1, POMT2, FKTN, ISPD, COL4A1, COL4A2, AP4M1, AP4S1, AP4E1, L1CAM, TREX1, RNASEH2A, RNASEH2B, RNASEH2C, RNASET2, SAMHD1, B3GALNT2, DYNC1H1, NSDHL, TMEM5, AP1S2, EML1

Sanger sequencing: ARX



!Extra MLPA analyse (SALSA P061 Lissencephaly probemix; MRC Holland)!



Targeted gene panel:

- ❖ Enrichment for a group of genes
- ❖ Beter coverage, cheaper, simple analysis (relative)
- ❖ No unsolicited findings
- ❖ Turn around time
- ❖ Disadvantage: self design, less flexibel than WES


Targeted capture/enrichment



- *Bait design: 103 genes > 1600 fragments*
- *(Agilent Sure Select*
- *Quality check (3 parameters)*
- *Sample prep/ Run Miseq; 12 patients/ run*
- *Mapping of data to the whole genome, average coverage ca. 300x*
- *Variant calling: SeqNext (JSI)*
- *Classification: Alamut*



!Extra: MLPA (SALSA P061 Lissencephaly; MRC Holland)



Classification sequence variants


- **Class 1; benign (frequent in control populations)**
 - Not reported
- **Class 2; silent variants and intronic variants, without effect on**
 - Not reported
- **Class 3; “Rest group (VUS)”**
 - Further investigation possible
- **Class 4; Probably pathogenic, HGMD link; earlier seen in patients**
 - Prenatal diagnosis and family testing (in agreement only)
- **Class 5; pathogenic (frameshift, nonsense, splice site (+/- 1 and 2))**
 - Prenatal diagnosis/ presymptomatic testing

Make a choice for which variant will be reported and/or confirmed by Sanger sequencing

Discussion **MEG**: Do genotype/ phenotype match (GENEMATCHER) → further investigation(s).

Wallis et al (ACGS & VKGL; 2013) Practice guidelines for the evaluation of pathogenicity and the reporting of sequence variants in Clinical Molecular Genetics


NGS-targeted: Pathogenic mutations (August 2013- Spring 2016)



Gene	Mutation	Inheritance
ACTB	c.625G>A, p.Val209Met	AD
ASPM	c.3202G>T, p.Glu1068*	AR
CENPJ	c.289dupA, p.Thr97fs	AR
CENPJ	c.1805_1808del, p.Glu602fs	AR
CEP152	c.4000_4001delAA, p.Lys1334fs	AR
CEP152	c.467dupG, p.Gln157fs	AR
COL18A1	c.3512_3515delCT, p.Leu1172Valfs	AR
COL18A1	c.2368C>T, p.Arg780*	AR
COL4A1	c.324+1G>A, p.?	AD
DCX	c.587G>A, p.Arg196His	X-Linked
DCX	c.907C>T, p.Arg303*	X-Linked
DYNC1H1	c.3359A>G, p.Tyr1120Cys	AD
DYNC1H1	c.3185A>C, p.Asp1062Ala	AD
ERCC2	c.1381C>G, p.Leu461Val	AR
FLNA	c.5160C>A, p.Tyr1720*	X-Linked
FLNA	c.4310_4314del, p.Pro1437fs	X-Linked
INTS8	c.1665_1699del5insT, p.Leu565Phefs	NN/AR
KIF2A	c.962A>C, p.His321Pro	AD
PAFAH1B1	c.1111C>T, p.Arg371*	AD
PAFAH1B1	c.337C>T, p.Arg113*	AD
PIK3CA	c.241G>A, p.Glu81Lys	AD
PIK3R2	c.1117G>A, p.Gly373Arg	AD
PNKP	c.1253_1269dup17, p.Thr424fs	AR
PNKP	c.1299-1G>T, p.?	AR
PNKP	c.793A>G, p.Met265Val	AR
PTEN	c.45A>C, p.Arg15Ser	AD
RAB3GAP1	c.748+1G>A, p.?	AR
RAB3GAP2	c.1361_1370del10, p.Asn454fs	AR
RAB3GAP2	c.265dupA, p.Met89fs	AR
RELN	c.4228G>A, p.Glu1410Lys	AR
RELN	c.6866C>T, p.Thr2289Ile	AR
RNASE2B	c.529G>A, p.Ala177Thr	AR
RNU4ATAC	g.122288508C>T, p.?	AR
RNU4ATAC	g.122288506C>T, p.?	AR
RNU4ATAC	g.122288445C>T, p.?	AR
RTTN	c.31+1G>T, p.?	AR
SLC25A9	c.720delC, p.Phe240fs	AR
TBC1D24	c.1008delT, p.His336fs	AR
TREX1	c.490dupA, p.Cys154fs	AD/AR
TSC2	c.4604A>C, p.Asp1535Ala	AD
TSEN54	c.3_8dup, p.Glu6_Pro7dup	AR
TSEN54	c.919G>T, p.Ala307Ser	AR
TUBA1A	c.791G>A, p.Arg264His	AD
TUBB2B	c.292G>A, p.Gly98Arg	AD
TUBB3	c.211G>A, p.Gly71Arg	AD
VPS13B	c.626-1G>C, p.?	AR

Total **46** pathogenic mutations; **22** likely pathogenic mutations in **32** different genes
> 600 VUS

Diagnostic yield: Targeted NGS





251 individuals with cortical brain malformations


Direct diagnosis: **18 (7,2%)**

After evaluation: **12 (4,8%)**

Total diagnostic yield: **12 %**

Diagnostic Yield depends on phenotyping






Why?

- ❖ *Patients genetically very heterogenous*
- ❖ *Same method for all panels: less money; less time*
- ❖ *Clinical Research Exome (CRE) (beter coverage disease genes)*
- ❖ *Diagnostic yield will increase*
- ❖ *Flexibel; genepanels can be expanded regularly (2x/ year)*
- ❖ *Open exome after (trio- analysis) is possible*
- ❖ *Whole Genome Sequencing Too expensive (at the moment)*

From 103 genes → 172 genes → **WES-CRE-NEUMIG**

Technical Information: WES-CRE NEUMIG




- ❖ Quality check
- ❖ *96 patients/1 run*
- ❖ Agilent SureSelect **Clinical Research Exome V2** capture (paired-end sequenced /Illumina platform (**outsourced**)).
- ❖ Average coverage of the exome is ~50x. Duplicate reads are excluded.
- ❖ Mapping to the genome: BWA-MEM algorithm (reference:<http://bio-bwa.sourceforge.net/>)
- ❖ Variant detection: Genome Analysis Toolkit HaplotypeCaller (reference:<http://www.broadinstitute.org/gatk/>)
- ❖ Filtering and annotation: Cartagenia software

- ❖ Classification with Alamut Visual



!Extra: MLPA (SALSA P061 Lissencephaly; MRC Holland)

WES-CRE: Pathogenic Mutations (Spring 2016)




Gene	Mutation	Inheritance
ANKLE2	c.1353+1G>C, p.?	AR
AP4M1	c.1317G>A, p.Trp439*	AR
AP4M1	c.1012C>T, p.Arg338*	AR
ASNS	c.1002T>G, p.Tyr334*	AR
ASPM	c.2495dupA, p.Tyr832*	AR
ASPM	c.8563A>T, p.Arg2855*	AR
CENPJ	c.3595C>T, p.Gln1199*	AR
CRB2	c.2887_2905del19, p.Met963fs	AR
ERCC2	c.1816G>T, p.Glu606*	AR
ERCC2	c.1381C>G, p.(Leu461Val)	AR
ERCC2	c.2150C>G, p.(Ala717Gly)	AR
ERMARD	c.1369_1373dup5, p.Leu459fs	AD
FLNA	c.6052G>T, p.Glu2018*	X-linked
ISPD	c.679C>T, p.Gln227*	AR
KIF11	c.1078_1085del, p.Leu3360Alafs*2	AD
KIF5C	c.2708G>A, p.(Arg903His)	AD
LAMA1	c.2755G>A, p.(Gly919Arg)	AR
LAMA1	c.2657C>T, p.(Ala886Val)	AR
LAMC	c.976+1G>A, p.?	AR
PAFAH1B1	c.430C>T, p.Arg144*	AD
PLK4	c.1299_1300ins20, p.Lys434*	AR
PLK4	c.1199_1200insGAA, p.His400delinsGlnAsp	AR
POMGNT1	c.880-1G>C, p.?	AR
POMGNT1	c.187C>T, p.Arg63*	AR
PTEN	c.355G>T, p.(Val119Phe)	AD
PTEN	c.518G>A, p.(Arg173His)	AD
PTEN	c.511C>T, p.Gln171*	AD
PTEN	c.959T>A, p.Leu320*	AD
RNASEH2B	c.529G>A, p.(Ala177Thr)	AR
RNASEH2B	c.136+1delG, p.?	AR
TREX1	c.309dupC, p.Thr104fs	AD
TSC2	c.1864C>T, p.(Arg622Trp)	AD
TUBA1A	c.641G>A, p.(Arg214His)	AD
TUBA1A	c.788C>T, p.(Pro263Leu)	AD
VLDLR	c.1962+2T>A, p.?	AR
VPS13B	c.8978A>G, p.(Asn2993Ser)	AR
WDR81	c.4668_4669delAG, p.Gly1557fs	AR

37 pathogenic / **8** likely pathogenic mutations in 25 different genes.

> **232** Variants of unknown significance

Diagnostic yield: WES-CRE




108 individuals with cortical brain malformations

Direct diagnosis: **8 (7,4%)**

After evaluation: **12(11,1%)**

Total diagnostic yield: **18,5 %**


Conclusions:



- Both approaches are suited (Advantages/ disadvantages)
- Feedback required for interpretation of VUS: genotype/ phenotype correlation
- Diagnostic yield improves after evaluation of results

Close collaboration of a Multidisciplinary Expert Group *is mandatory for the diagnosis of complex MCD patients*

Acknowledgements:



Department Clinical Genetics

Diagnostic Laboratory

Marjon van Slegtenhorst
Marianne van Tienhoven



Galhana Bolman
Joan Kromosoeto

Rick van Minkelen

Walter de Valk

Dicky Halley

Counseling:

Grazia Mancini

Alice Brooks
(Maartje Berghuis)
(Renske Oegema)

Rachel Schot

Functional unit:

Frans Verheijen

Medical expertise team:

MEG-NEUMIG



Sequencing:

Dept. of Biomics, Erasmusmus MC, The Netherlands
GenomeScan B.V., Leiden, The Netherlands



Q&A